Product Specification

Disinfecting Wipes Brand: Mediwipes





If you would like samples or have any questions about our Mediwipes, please email michelle@finetouch.co.nz or call us on 0800 82 82 82

Product Features

Single use disinfecting wipe

Wipe size: 200 x 140mm

Textured wipe - helps remove dirt, blood, organic matter and salvia from surfaces

Perforated wipe – ensures no wastage when dispensed from container

Pop up dispenser shrink wrapped – helps retain moisture

Moist wipe - ensures surface is sufficiently covered by the formulation for product efficacy

160 wipes per dispenser, 12 dispensers per carton

Product Code: C-MEDIW (carton)

Product Ingredients

Alcohol 70%

Chlorhexidine digluconate 0.50%

Centrimonium bromide 0.16%

Water Balance

Product Specification

Disinfecting Wipes Brand: Mediwipes



Quality Assurance

ISO 9001: 2008 certified by SGS (UK)

NZ Medical Device Class 1: WAND Reference 091029 - WAND - 67R7VW

MPI Approved C44 for food service

MPI Approved for use in farm dairies

Asure Quality Assessed for food, beverage, dairy factory, food contact surfaces (ref H2428)

Approved cleaning product by Howard Wright for medical beds and stretchers, mattress covers, upholstery tops and shower trolley tub fabrics

99% bactericidal effect – Results from NZ Laboratory Services Ltd on the following organisms tested for a 2 minute kill time at 100%

Staphylococcus aureus

Escherichia coli

Listeria monocytogenes

Salmonella Menston

Pseudomonas aeruginosa

Methicillin Resistant Staphylococcus aureus (MRSA)

Streptococcus dysgalactiae

Streptococcus agalactiae

Streptococcus uberis

Product Uses

Mediwipes are ideal for:

Hospitals

Medical clinics

Dental clinics

Laboratories

Vet clinics

Offices

Schools

Food processing areas

Food service industry (restaurant, cafe)

Food probe disinfection

Fine Touch Disposables, PO Box 11 312, Christchurch 8443. Ph 03 348 8516, 0800 82 82 82.

www.finetouch.co.nz